DATASHEET

The Adult Deluxe Airway Management Trainer allows training of all standard techniques for airway maintenance with the possibility of laryngospasm and tongue edema as well as training of cricothyrotomy.

All anatomical landmarks are recognizable and the vocal cords are in a neutral position and highlighted for easy viewing. The lifelike skinis very stretchable and robust, and the upper teeth break out if inappropriate technique is applied during intubations. To determine a successful intubation both lungs and stomach are visible.

The cricoid area has replaceable skins for an economical way to teach and practice surgical airway techniques.

The trainer is latex free and without carotid pulse.

The board is anchored to surfaces by suction cup feet.

This head is also available without the board to use as a replacement for the SMART STAT head (101-459U Head for PDA STAT Simulation, with carotid pulse).

Adult Deluxe Airway Management Trainer with Board



KEY BENEFITS

- Cricothyrotomy
- Laryngospasm and tongue edema
- All anatomical landmarks
- Vocal cords in neutral position
- Upper teeth break out
- Lifelike skin and flexible tongue

ORDERING SPECIFICATIONS

Item no.	
101-502U	Adult Deluxe Airway Management Trainer with Board
101-459U	Head for PDA STAT Simulation



■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com