

DATASHEET

This Obstetrical Manikin is a lifelike birth simulator including an anatomically correct pelvic model with newborn and placenta allowing your trainees to practice normal vaginal birth as well as procedures for emergency childbirth.

The simulator contains two interchangeable abdominal covers - one is a modular pregnant belly overlay with permanently installed fetus for training practice of Leopold's maneuver to determine the fetus' lie by palpating its skull and kneecaps, and the other is a clear abdominal overlay to see position during labor.

The pelvic cavity has pronounced pelvic landmarks, located spinal column, angled birth canal, ilium, ischium, sacrum, sacro spinous ligaments, and greater sciatic notch.

The fetus provided is of a realistic newborn size with fontanels and cranial sutures.

Obstetrical Manikin



KEY BENEFITS

- Lifelike birth simulator with newborn and placenta
- Realistic pelvic cavity
- Interchangeable abdominal covers
- Training of Leopold's maneuver
- Fetus with fontanels and cranial sutures

ORDERING SPECIFICATIONS

Item no.

Obstetrical Manikin

110-180U
incl. disposable umbilical cords with clamps, extra vulva, blood powder, modular pregnant belly overlay with permanently installed fetus, clear abdominal overlay, soft carry bag

110-182U
Obstetrical vulva

110-185U
Placenta

110-159U
Umbilical cord

Ambu

■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com