

STAT Baby

DATASHEET

The STAT Baby is a lifelike simulator with multiple features for PALS training and represents a 9-month-old infant with lifelike skin and a weight of 8 kg.

It features advanced airway management, I.V. access, PICC access, I.O. access, chest decompression, chest tube placement, as well as urinary catheterization, rectal medication, ECG monitoring and defibrillation, intramuscular and subcutaneous injections, and more.

Airway management

This simulator has an advanced airway allowing the practice of ET and laryngeal tubes, laryngeal masks, oropharyngeal and nasopharyngeal airways. The advanced intubation head enables also tongue edema and laryngospasm as well as the Sellick's maneuver. In addition, the placement of a nasogastric tube for gavage or lavage can be trained.

Pulses

It features seven manual pulse points - pairs at carotid, femoral and brachial, plus the left radial pulse.

Intravenous access

The intravenous access is possible on the right side antecubital, left dorsal hand and left saphenous pedal. PICC line access is possible on the right antecubital fosse.

Intraosseous leg

The baby allows practice of intraosseous infusion on the right leg.

Chest tube placement & chest decompression

The baby features also bilateral tension pneumothorax reduction, and chest tube placement with fluid return on the left midaxillary site. PICC line access is possible on the right antecubital fosse.

Interactive ECG Simulator

With the Interactive ECG Simulator 17 pediatric heart rhythms can be recreated, as well as simulated pacing or cardioversion. It can be used with manual, semi-automatic or automatic defibrillators.



KEY BENEFITS

- Advanced intubation head with nasal intubation, tongue edema and laryngospasm
- 4-lead monitoring and defibrillation
- I.O. access, right leg
- Peripheral I.V. access, 3 sites on hand and foot
- PICC access site
- Intramuscular and subcutaneous injections
- 7 pulses, manual
- Nasogastric tube placement
- Tension pneumothorax reduction
- Chest tube placement
- Urinary catheterization, female
- Rectal medication

ORDERING SPECIFICATIONS

Item no.

101-350U	STAT Baby Fullbody infant manikin with Interactive ECG Simulator, defibrillator adapters, IV vein replacement kit (1 arm and 1 leg), replacement arm and leg skins, 12 IO leg bones, 4 pneumothorax pads, 2 IM injection pads, 4 reservoir bags, blood powder, pulse bulb, syringes, carry bag
101-095U	Arm Vein Replacement Kit

Ambu

■ Ambu A/S
Baltorpbakken 13
2750 Ballerup
Denmark
T +45 72 25 20 00
F +45 72 25 20 50
ambu.com